

Single-Mode Digital U-Type Module



IPS proprietary Wavelength Stabilized Laser features high output power with narrow spectral bandwidth. The laser's stabilized peak wavelength remains "locked" regardless of case temperature (15 to 45 °C). Devices can be spectrally tailored to suit application needs and offer side mode suppression ratios (SMSRs) better than 40 dB, thereby providing extremely high signal to noise ratio and making these sources ideal for Raman spectroscopy and pump laser applications. The laser is integrated with high performance laser drive and temperature control electronics in a compact package. In addition to integration into systems, IPS' OEM U-Type module is designed to "drop in" to our UL/CE and IEC certified turnkey modules to offer wavelength flexibility at a lower cost.

Applications

This laser package is designed for OEM Integration and is ideal for:

- High Resolution Raman Spectroscopy
 - Confocal Microscopy
 - Raman Imaging
 - Portable Raman
 - Process Raman
- Direct-diode Frequency Doubling
- Fiber Laser Seeding
- Metrology & Interferometry
- Remote Sensing

Key Features

- High-Power Single-Spatial-Mode, Single-Frequency Output
- Narrow Spectral Linewidth (< 100 MHz FWHM)
- High Power Single-mode Fiber Coupled Output
- Excellent Beam Quality ($M^2 < 1.1$)
- Integral ESD Protection & Thermistor
- Temperature Stabilized Spectrum (< 0.007 nm/°C)
- Low Power consumption (< 5.5 W)
- > 45 dB SMSR Typical
- 3" x 2.5" x 0.69" Package Weighing < 4 oz

Standard Wavelengths

633nm	785nm	1030nm
638nm	808nm	1053nm
780nm	830nm	1064 nm
783nm	976nm	


Specifications



Wavelegnth Tolerance	+/- 0.5nm
Spectral Linewidth	<100MHZ
Wavelength Stability Range	15 - 45 °C
SMSR	45 - 55 dB typical
Polarization Extinction (PER)	>17 dB (20 dB Typical)
Polarization Orientation	Standard is PM slow
Output Power Stability	1% typical
Modulation Rate	CW to 1KHz (for 10% power to CW) up to 10 KHz for 50% power.
Warm-Up Time	10 seconds from cold start
	1.5 seconds from warm start

λ (nm)	Min. Power (mW)	Base Part Number
633	20	I0633SU0020PA-USB
638	25	I0638SU0025PA-USB
780	50	I0780SU0050PA-USB
783	50	I0783SU0050PA-USB
785	50	I0785SU0050PA-USB
	100	I0785SU0100PA-USB
808	50	I0808SU0050PA-USB
830	50	I0830SU0050PA-USB
976	220	I0976SU0220PA-USB
1053	50	I1053SU0050PA-IS-USB
	120	I1053SU0120PA-USB
	300	I1053SU0300PA-USB
1064	50	I1064SU0050PA-IS-USB
	120	I1064SU0120PA-USB
	300	I1064SU0300PA-USB

Part Schema




Selected Data



Supply Voltage	4.9V min to 5.1V max
Power Consumption	3.5W typical, 5.5W maximum
Photodiode Current	30 uA
Laser Set Point Control (LD SET)	900mA to 1000mA when pin 2 grounded

Optical Fiber Options	Single-Mode Fiber
	Polarization Maintaining, Panda Type
Connector	FC/APC
Electrical Connector	10-pin, Molex #53014-1010 (mating connector: 51004-1000)
Module Dimensions	3.0 x 2.5 x 0.69 inches
Module Weight	100 grams (3.5 ounces)
Case Material	Anodized Aluminum
Operating Temperature	10 to 45 °C
Cooling air flow (internal)	100 LRM with attached heatsink
Environment	0-80% Humidity, non-condensing
Storage Temperature	-10 to 55 °C



Custom Capability

- Custom wavelengths available upon request
- FC/PC, FC/APC, or SMA output coupler
- Single-mode or Polarization-maintaining fiber available with orientation in either fast or slow axis
- Integral optical isolator available
- External TEC (e.g. No TEC inside of package optional)

Electrical Specs





Pin	Pin Label	Function
1*	V+	Power Supply: 5V-12V. For lasers operating <600mW, a 5V minimum is required. For lasers operating >600mW, a minimum of 9V is required, 12V recommended
2^	VBIAS ENABLE (Set Enable)	For analog control: Enable external laser power control through pin 8 (LD VBIAS): High (>3.6V up to V+ supply voltage) = Enable Low (GND) = Disable The same Functionality is emulated via software as a replacement for this analog control.
3	SCL/TX	I2C: SCL standard (RS232: TX; board by request)
4	SDA/RX	I2C: SDA standard (RS232: RX; board by request)
5**	GND	Ground
6*	V+	See Pin 1
7	Enable	Laser Enable: TTL High (>3.6V up to V+ supply voltage) = Enable TTL Low (GND) = Disable Pin may also be used for PW modulation at a rate of 1kHz.
8^	LD VBIAS (LD Set)	Laser power setpoint - Enables analog external control of laser drive current. Drive VBIAS Enable (pin2) high to enable this option. Drive voltage between 0V and 1V. The voltage bias will be a 1:1 ratio to the laser drive current. See Operational Note #6. ..
9^	PD+	For analog readout. Connect voltmeter to PD+ (pin 9) and GND for photo diode V output (0V-3.3V). The same functionality is emulated via software as a replacement for this analog control.
10**	GND	Ground

Pins 1, 2, 6, 7 and 10 are required for laser operation

Notes: *Power must be supplied to both V+ pins (pin 1 and pin 6), **GND must be supplied to both GND pins

(pin 5 and pin 10). ^ Pins 2, 8, and 9 are optional for analog control/readout

Mechanical Drawings



- Connector Options:
1. FC/PC (M8 Thread)
 2. SMA 905 (1/4"-36 Thread)

Operational Notes

1. To adjust power output, IPS recommends Pulse Width Modulation (PWM) to adjust AVERAGE power rather than using pin 8 (LD SET) for single-mode diode lasers. See Note 2.
2. By using PWM, user can adjust average power from 10% to 100% in digital increments by setting pulse width and duty cycle. For example, if a 50% duty cycle is selected, the laser will be on 50% of the time, and off 50% of the time, making the average power equal to 50% of the CW output power. The sample will experience a lower average power. Rise/fall time is approximately 20 microseconds.
3. Heat sink and 5V power supply are not included with module. Please ask about our turn-key package that is available as an add-on.
4. Do not retro-reflect beam! This can cause Catastrophic Optical Damage (COD) and is not covered under warranty
5. Be aware that this approach may cause laser mode hopping behavior in single-mode lasers. Do not exceed maximum recommended drive current for given laser model. The same functionality is emulated via software as a replacement for this analog control.

Innovative Photonic Solutions, Inc.
313 Enterprise Drive
Plainsboro, NJ 08536

Phone: (732) 355-9300

sales@ipslasers.com
www.ipslasers.com

